

X-Ray Fluorescence

March 13, 2018

1 X-Ray Fluorescence

In this lab, we will use an Amptek Mini-X X-Ray Tube in conjunction with the X-123 Complete Spectrometer to use X-Ray fluorescence to determine the composition of an unidentified sample.

1.1 6 March

We were unable to do any work on the experiment last week, as we were trying to modify the fluorescence setup to do X-Ray diffraction. However, after consulting with Jun Ye, we decided to abandon the diffraction experiment in favor of performing the fluorescence experiment.

SETUP * We first had to reinstall the device drivers for the spectrometer. All driver files downloaded to desktop. * Once the computer is talking to the devices, we were able to make a run on a Pb sample. * We follow the procedures given by the Experimenter's Kit Quick Start June 2012.pdf

CALIBRATION

- We have several different labeled samples to calibrate our equipment with: Cu, Pb, Sn, Al and we assume Steel.
- According to the X-123 Specification sheet, the Energy Range at $> 25\%$ is 1.5 keV to 25 keV, although the manual claims it will detect outside this range at lower energy (<http://amptek.com/products/x-123-complete-x-ray-spectrometer-with-si-pin-detector/>)
- Combining this range with the K and L emission lines (<http://amptek.com/pdf/xraychrt.pdf>), we can use our samples to calibrate the energy axis of the detector. We will do this for each of the samples to get an idea of the precision and accuracy of the spectrometer.
- We follow the following calibration procedure:
 - Place Sample in spectrometer
 - Clear all data from Amptek DppMCA (Gain set to 30), time duration set to 30 s
 - Set X ray tube to 30 kV, 100 μA
 - Start X ray
 - Clear Spectrometer Data
 - Start Spectrometer collection
 - Turn off xray when data collection complete.
 - Compare sample peaks to range of spectrometer to determine which peaks correspond to which transistions

- Calibrate Spectrometer
- Note: To use the calibration feature of the Amptek software, we first define Regions of Interest (ROI) by providing the min/max channel values of the windows defining the peaks.
- We then go into calibration, click on the peak of interest, click 'centroid' on the RHS, and enter the actual value in keV of the peak location. Then click 'add'
- Repeat for all peaks in the spectrum
- Once complete, click 'enable calibration.'

We repeat this procedure 3 times for each of our samples.

* A = -0.0311604 keV (Offset) * B = 0.0176715 keV (keV/Channel)

1.2 7 March

We had difficulty with the Analysis software from Amptek. We will repeat the calibration, and use the DppMCA interface with the spectrometer for data collection.

- We set the Voltage and current of the X-Ray tube to 40 kV and 100 μ A
- We collect the spectra of Pb, Cu, and Sn for calibration.
- We calibrate the horizontal scale based on the known energies of the various known transition lines.

Through calibration, we have a way of characterizing the error in our spectrometer, as well as converting an arbitrary channel into an energy. Once this is complete, we have a way of determining the energies of peaks in various samples and comparing them to the 'fingerprints' of the constituent atoms. This should allow us to determine the composition of various samples.

```
In [31]: import pandas as pd
import numpy as np
import matplotlib.pyplot as plt
import os
from scipy import signal
from scipy import stats
from scipy import optimize
from ipywidgets import interact
from scipy.stats import norm

plt.style.use('ggplot')
plt.rc('text', usetex=True)
plt.rc('font', family='serif')
plt.rcParams.update({'font.size': 18, 'text.usetex': True})

%matplotlib inline

In [61]: class spectrograph:
        """
```

From a datafile containing the counts from an X-Ray Fluorescence Spectrograph, we read the counts from a spectrometer for energies in a given channel. See the 'MCAC' line for the specific number of channels used. The ones we used here were for 8192 channels

```

def __init__(self, file_path):
    self.file_path = file_path
    self.file = pd.read_csv(file_path, skiprows = 32, skipfooter = 71, names = ['Co
    self.cal_data = pd.read_csv(self.file_path, sep = ' ', skiprows = 13, skipfooter
                                names = ['Channel', 'Energy (keV)'], engine = 'pytho
    self.cal_fit_params, cal_fit_covar = np.polyfit(self.cal_data['Channel'], self.
                                                1, cov=True)

    self.cal_err = np.sqrt(np.diag(cal_fit_covar))
    self.calm_report = self.reported_values(self.cal_fit_params[0], self.cal_err[0])
    self.calb_report = self.reported_values(self.cal_fit_params[1], self.cal_err[1])
    self.cal_energy = np.poly1d(self.cal_fit_params)(np.arange(len(self.file)))
    self.cal_chi2, self.cal_p = stats.chisquare(self.cal_data['Energy (keV)'], \
                                                np.poly1d(self.cal_fit_params)(self
    self.file.insert(0, 'Energy (keV)', self.cal_energy)

def print_data(self):
    return self.file

def get_counts(self):
    return self.file['Counts']

def reported_values(self, value, error):
    a = int(np.floor(np.log10(np.abs(error))))
    report_Err = round(error, -a)
    report_Val = round(value, -a)
    return report_Val, report_Err

def plot_calib(self):
    cal_fit_x = np.linspace(0,9000,9000)
    cal_fit_y = np.poly1d(self.cal_fit_params)(cal_fit_x)
    fig, ax = plt.subplots(figsize = (10,10))
    ax.scatter(self.cal_data['Channel'], self.cal_data['Energy (keV)'], color = 'b')
    ax.plot(cal_fit_x, cal_fit_y, \
            label = str('Best Fit:\n m = {}' + r'$\pm$' + '{} keV/Channel \n b = {}'
                        + r'$\pm$' + '{} keV\n p-val = {}')) \
            .format(self.calm_report[0], self.calm_report[1], \
                    self.calb_report[0], self.calb_report[1], self.cal_p))
    ax.set_title('Channel to Energy Calibration')
    ax.set_xlabel('Channel Number')
    ax.set_ylabel('Energy (keV)')
    ax.axhline(y=0, color = 'k')
    ax.axvline(x=0, color = 'k')

```

```

        ax.legend()
        fig.tight_layout()

    def plot_spectrum(self, dataName):
        self.fig, self.ax = plt.subplots(figsize = (10,10))
        self.ax.plot(self.file['Counts'])
        self.ax.plot(self.file['Energy (keV)', self.file['Counts'])
        self.ax.set_title('Channel Counts for {}'.format(str(dataName)))
#         self.ax.set_xlabel('Channel (a.u)')
        self.ax.set_xlabel('Energy (keV)')
        self.ax.set_ylabel('Counts (a.u.)')
        self.ax.axhline(y=0, color = 'k')
        self.ax.axvline(x=0, color = 'k')
        self.fig.tight_layout()

    def find_peaks(self, peak_width, snr):
        self.peaks = self.file['Energy (keV)'][signal.find_peaks_cwt(self.file['Counts'],
                                                                    np.arange(1,peak_w

    def fit_peaks(self, peaklocs, peakhts, peakwd):
        peak_info = pd.DataFrame(np.zeros(len(peaklocs),6),\
                                columns = ['Peak Location', 'Peak Width', 'Peak Height',\
                                             'Peak Location Err', 'Peak Width Err', 'Peak Height Err'])
        for i in np.arange(len(peaklocs)):
            optimize.curve_fit(gauss_fit, )


    def gauss_fit(x, x0, height, width):
        return height*width*np.sqrt(2*np.pi)*norm.pdf(x, loc = x0, scale = width)

    def make_peak(x0, ):
        x = np.linspace(0, 9000,9000)
        y = height*np.sqrt(2*np.pi)*norm.pdf(x, loc = x0, scale = )
        fig, ax = plt.subplots(figsize = (8,8))
        ax.plot(x,y)
        fig.tight_layout()


In [33]: a = spectrograph('SumCrossXRF/Calibration/7Mar_Cu1.txt')
a.find_peaks(100, 20)
a.plot_spectrum('Copper Calibration, Run 1')
a.ax.vlines(a.peaks, 0, 3500, color = 'b')
# a.peaks
# np.poly1d(a.cal_fit_params)(100)

Out[33]: <matplotlib.collections.LineCollection at 0x7f16d67952b0>

```





```
In [46]: a.plot_calib()
```





```
In [64]: test = spectrograph('SumCrossXRF/Calibration/7Mar_Cu1.txt')
# test.cal_good_fit()
# print(test.calm_report[1])
# test.plot_calib()
# test.file
# len(test.cal_energy)
# test.file['Energy (keV)']
# np.array(test.file['Counts'])


In [57]: # Plot all calibration Runs
for file in sorted(os.listdir('SumCrossXRF/Calibration/')):
    spect = spectrograph('SumCrossXRF/Calibration/{0}'.format(file))
    spect.plot_spectrum(file[5:7])
```



















```
In [62]: # Plot all Data Runs
for file in sorted(os.listdir('SumCrossXRF/Data/')):
    #     print(file.partition('.')[0].partition('_')[2])
    spect = spectrograph('SumCrossXRF/Data/{}'.format(file))
    spect.plot_spectrum(file.partition('.')[0].partition('_')[2])
```




































```
In [37]: # Plot all copper data runs for width comparisons
for file in ['SumCrossXRF/Data/7Mar_0-175mmCu.txt', 'SumCrossXRF/Data/7Mar_0-20mmCu.txt']:
    spect = spectrograph(file)
    spect.plot_spectrum(file.partition('.')[0].partition('_')[2])
```







Goal: Do a fitting: If an element shows up in a material, *all* peaks should be showing up in the material spectrum. Create the theoretical spectrum for each element by creating gaussians centered on the theoretical positions. First, we parse data from <http://www.med.harvard.edu/jpnm/physics/refs/xrayemis.html> to get the location of all x-ray fluorescence peaks. Then, we will use an arbitrary data run from the instrument to generate a plot of these peaks against the same bins as the instrument collected data. We assume a gaussian fit for each peak.

```
In [38]: wavelengths = pd.read_fwf('X-Ray_database', names = ['Z', 'Ka1', 'Ka2', 'Kb1', 'La1', 'Lg1'], index_col = 0, skiprows = 2, usecols = [0,2,3,4,5,6,7,8,9]).T
wavelengths=wavelengths.fillna(0.)
wavelengths=wavelengths.replace(['-', '0'],0.)
wavelengths
for i in wavelengths:
```

```
wavelengths[i] = wavelengths[i].astype(float)
wavelengths.columns.astype('int')
# len(wavelengths.columns)
wavelengths
```

```
Out[38]: Z      3      4      5      6      7      8      9     10     11  \
Ka1  0.0543  0.1085  0.1833  0.277  0.3924  0.5249  0.6768  0.8486  1.04098
Ka2  0.0000  0.0000  0.0000  0.000  0.0000  0.0000  0.0000  0.8486  1.04098
Kb1  0.0000  0.0000  0.0000  0.000  0.0000  0.0000  0.0000  0.0000  1.07110
La1  0.0000  0.0000  0.0000  0.000  0.0000  0.0000  0.0000  0.0000  0.00000
La2  0.0000  0.0000  0.0000  0.000  0.0000  0.0000  0.0000  0.0000  0.00000
Lb1  0.0000  0.0000  0.0000  0.000  0.0000  0.0000  0.0000  0.0000  0.00000
Lb2  0.0000  0.0000  0.0000  0.000  0.0000  0.0000  0.0000  0.0000  0.00000
Lg1  0.0000  0.0000  0.0000  0.000  0.0000  0.0000  0.0000  0.0000  0.00000

Z      12  ...      86      87      88      89      90  \
Ka1  1.2536  ...   83.7800  86.1000  88.4700  90.8840  93.3500
Ka2  1.2536  ...   81.0700  83.2300  85.4300  87.6700  89.9530
Kb1  1.3022  ...   94.8700  97.4700  100.1300  102.8500  105.6090
La1  0.0000  ...   11.7270  12.0313  12.3397  12.6520  12.9687
La2  0.0000  ...   11.5979  11.8950  12.1962  12.5008  12.8096
Lb1  0.0000  ...   14.3160  14.7700  15.2358  15.7130  16.2022
Lb2  0.0000  ...    0.0000  14.4500  14.8414    0.0000  15.6237
Lg1  0.0000  ...   16.7700  17.3030  17.8490  18.4080  18.9825

Z      91      92      93      94      95
Ka1  95.8680  98.4390  0.0000  0.0000  0.0000
Ka2  92.2870  94.6650  0.0000  0.0000  0.0000
Kb1  108.4270  111.3000  0.0000  0.0000  0.0000
La1  13.2907  13.6147  13.9441  14.2786  14.6172
La2  13.1222  13.4388  13.7597  14.0842  14.4119
Lb1  16.7020  17.2200  17.7502  18.2937  18.8520
Lb2  16.0240  16.4283  16.8400  17.2553  17.6765
Lg1  19.5680  20.1671  20.7848  21.4173  22.0652
```

[8 rows x 93 columns]

```
In [63]: spectral_data[spectral_data.columns[1:]];
```

```
In [48]: @interact(x0 = (0,9000,50),  = (0,1000,20),  continuous_update=False)
```

```
def make_peak(x0, ):
    x = np.linspace(0, 9000,9000)
    y = *np.sqrt(2*np.pi)*norm.pdf(x, loc = x0, scale = )
    fig, ax = plt.subplots(figsize = (8,8))
    ax.plot(x,y)
    fig.tight_layout()
```

```
interactive(children=(IntSlider(value=4500, description='x0', max=9000, step=50), IntSlider(valu
```

```

In [41]: # Create a dataset of the spectrum discretized into the same size as the instrument col
def make_spectrum(x, wavelengths, ):
    spectrum = np.zeros(len(x))
    for x0 in wavelengths:
        add = *np.sqrt(2*np.pi)*norm.pdf(x, loc = x0, scale = )
        if (x0 <1.5 or x0>25.):
            spectrum += .01*np.sqrt(2*np.pi)*norm.pdf(x, loc = x0, scale = )
        else:
            spectrum += *np.sqrt(2*np.pi)*norm.pdf(x, loc = x0, scale = )
    return spectrum

# x = np.linspace(0,40,9000)
# # .5*np.sqrt(2*np.pi)*norm.pdf(x, loc = )
# fig, ax = plt.subplots(figsize = (10,10))
# ax.plot(x, make_spectrum(x, wavelengths[87], .1))

spectral_data = {}
for column in wavelengths:
    spectral_data[column] = make_spectrum(test.file['Energy (keV)'],wavelengths[column])

spectral_data = pd.DataFrame(spectral_data)
spectral_data.insert(0,'Energy (keV)',test.file['Energy (keV)'])
spectral_data
# len(spectral_data)

# fig, ax = plt.subplots()
# ax.plot(spectral_data['Energy (keV)'],spectral_data[29])
# spectral_data[5]

```

```

Out[41]:

```

| | Energy (keV) | 3 | 4 | 5 | 6 | 7 \ |
|----|--------------|----------|----------|----------|----------|----------|
| 0 | 0.045407 | 0.073104 | 0.071338 | 0.067008 | 0.063828 | 0.063167 |
| 1 | 0.049811 | 0.071823 | 0.070251 | 0.065936 | 0.062590 | 0.061861 |
| 2 | 0.054214 | 0.070433 | 0.069063 | 0.064780 | 0.061269 | 0.060466 |
| 3 | 0.058617 | 0.068941 | 0.067781 | 0.063547 | 0.059872 | 0.058989 |
| 4 | 0.063020 | 0.067355 | 0.066410 | 0.062244 | 0.058406 | 0.057437 |
| 5 | 0.067424 | 0.065682 | 0.064959 | 0.060878 | 0.056880 | 0.055819 |
| 6 | 0.071827 | 0.063932 | 0.063434 | 0.059456 | 0.055303 | 0.054143 |
| 7 | 0.076230 | 0.062112 | 0.061842 | 0.057987 | 0.053682 | 0.052417 |
| 8 | 0.080634 | 0.060232 | 0.060192 | 0.056476 | 0.052027 | 0.050650 |
| 9 | 0.085037 | 0.058300 | 0.058489 | 0.054932 | 0.050345 | 0.048850 |
| 10 | 0.089440 | 0.056325 | 0.056743 | 0.053361 | 0.048646 | 0.047025 |
| 11 | 0.093844 | 0.054316 | 0.054961 | 0.051770 | 0.046936 | 0.045184 |
| 12 | 0.098247 | 0.052281 | 0.053149 | 0.050166 | 0.045225 | 0.043334 |
| 13 | 0.102650 | 0.050229 | 0.051315 | 0.048556 | 0.043520 | 0.041482 |
| 14 | 0.107054 | 0.048168 | 0.049466 | 0.046945 | 0.041827 | 0.039638 |
| 15 | 0.111457 | 0.046106 | 0.047609 | 0.045339 | 0.040154 | 0.037807 |
| 16 | 0.115860 | 0.044051 | 0.045750 | 0.043743 | 0.038507 | 0.035996 |
| 17 | 0.120264 | 0.042010 | 0.043896 | 0.042163 | 0.036893 | 0.034211 |

| | | | | | | |
|------|-----------|----------|----------|----------|----------|----------|
| 18 | 0.124667 | 0.039989 | 0.042052 | 0.040603 | 0.035316 | 0.032460 |
| 19 | 0.129070 | 0.037995 | 0.040224 | 0.039066 | 0.033782 | 0.030745 |
| 20 | 0.133473 | 0.036034 | 0.038417 | 0.037557 | 0.032294 | 0.029074 |
| 21 | 0.137877 | 0.034110 | 0.036636 | 0.036078 | 0.030858 | 0.027450 |
| 22 | 0.142280 | 0.032231 | 0.034885 | 0.034633 | 0.029475 | 0.025878 |
| 23 | 0.146683 | 0.030398 | 0.033169 | 0.033223 | 0.028150 | 0.024360 |
| 24 | 0.151087 | 0.028617 | 0.031490 | 0.031851 | 0.026883 | 0.022901 |
| 25 | 0.155490 | 0.026891 | 0.029852 | 0.030518 | 0.025677 | 0.021502 |
| 26 | 0.159893 | 0.025222 | 0.028259 | 0.029226 | 0.024533 | 0.020166 |
| 27 | 0.164297 | 0.023614 | 0.026711 | 0.027974 | 0.023452 | 0.018894 |
| 28 | 0.168700 | 0.022067 | 0.025212 | 0.026764 | 0.022433 | 0.017689 |
| 29 | 0.173103 | 0.020584 | 0.023763 | 0.025595 | 0.021476 | 0.016550 |
| ... | ... | ... | ... | ... | ... | ... |
| 8161 | 35.980840 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8162 | 35.985243 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8163 | 35.989646 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8164 | 35.994050 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8165 | 35.998453 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8166 | 36.002856 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8167 | 36.007260 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8168 | 36.011663 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8169 | 36.016066 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8170 | 36.020469 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8171 | 36.024873 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8172 | 36.029276 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8173 | 36.033679 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8174 | 36.038083 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8175 | 36.042486 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8176 | 36.046889 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8177 | 36.051293 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8178 | 36.055696 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8179 | 36.060099 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8180 | 36.064503 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8181 | 36.068906 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8182 | 36.073309 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8183 | 36.077713 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8184 | 36.082116 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8185 | 36.086519 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8186 | 36.090922 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8187 | 36.095326 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8188 | 36.099729 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8189 | 36.104132 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |
| 8190 | 36.108536 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | 0.000000 |

| | 8 | 9 | 10 | 11 | ... | 86 | 87 | 88 | \ |
|---|----------|----------|----------|----------|-----|----------|----|----|---|
| 0 | 0.063143 | 0.063143 | 0.054123 | 0.045102 | ... | 0.009020 | 0 | 0 | |
| 1 | 0.061833 | 0.061833 | 0.053000 | 0.044167 | ... | 0.008833 | 0 | 0 | |
| 2 | 0.060433 | 0.060433 | 0.051800 | 0.043167 | ... | 0.008633 | 0 | 0 | |

| | | | | | | | | |
|------|----------|----------|----------|----------|-----|----------|----|----|
| 3 | 0.058951 | 0.058950 | 0.050529 | 0.042107 | ... | 0.008421 | 0 | 0 |
| 4 | 0.057393 | 0.057393 | 0.049194 | 0.040995 | ... | 0.008199 | 0 | 0 |
| 5 | 0.055768 | 0.055768 | 0.047801 | 0.039834 | ... | 0.007967 | 0 | 0 |
| 6 | 0.054084 | 0.054084 | 0.046358 | 0.038631 | ... | 0.007726 | 0 | 0 |
| 7 | 0.052350 | 0.052349 | 0.044871 | 0.037392 | ... | 0.007478 | 0 | 0 |
| 8 | 0.050573 | 0.050572 | 0.043348 | 0.036123 | ... | 0.007225 | 0 | 0 |
| 9 | 0.048762 | 0.048761 | 0.041795 | 0.034829 | ... | 0.006966 | 0 | 0 |
| 10 | 0.046924 | 0.046923 | 0.040220 | 0.033517 | ... | 0.006703 | 0 | 0 |
| 11 | 0.045069 | 0.045068 | 0.038629 | 0.032191 | ... | 0.006438 | 0 | 0 |
| 12 | 0.043202 | 0.043201 | 0.037030 | 0.030858 | ... | 0.006172 | 0 | 0 |
| 13 | 0.041334 | 0.041332 | 0.035428 | 0.029523 | ... | 0.005905 | 0 | 0 |
| 14 | 0.039469 | 0.039467 | 0.033829 | 0.028191 | ... | 0.005638 | 0 | 0 |
| 15 | 0.037615 | 0.037614 | 0.032240 | 0.026867 | ... | 0.005373 | 0 | 0 |
| 16 | 0.035780 | 0.035777 | 0.030666 | 0.025555 | ... | 0.005111 | 0 | 0 |
| 17 | 0.033968 | 0.033965 | 0.029113 | 0.024261 | ... | 0.004852 | 0 | 0 |
| 18 | 0.032185 | 0.032182 | 0.027584 | 0.022987 | ... | 0.004597 | 0 | 0 |
| 19 | 0.030437 | 0.030433 | 0.026086 | 0.021738 | ... | 0.004348 | 0 | 0 |
| 20 | 0.028729 | 0.028724 | 0.024621 | 0.020517 | ... | 0.004103 | 0 | 0 |
| 21 | 0.027064 | 0.027058 | 0.023193 | 0.019327 | ... | 0.003865 | 0 | 0 |
| 22 | 0.025446 | 0.025440 | 0.021806 | 0.018171 | ... | 0.003634 | 0 | 0 |
| 23 | 0.023880 | 0.023872 | 0.020461 | 0.017051 | ... | 0.003410 | 0 | 0 |
| 24 | 0.022366 | 0.022357 | 0.019163 | 0.015969 | ... | 0.003194 | 0 | 0 |
| 25 | 0.020909 | 0.020898 | 0.017912 | 0.014927 | ... | 0.002985 | 0 | 0 |
| 26 | 0.019509 | 0.019496 | 0.016711 | 0.013926 | ... | 0.002785 | 0 | 0 |
| 27 | 0.018168 | 0.018153 | 0.015560 | 0.012966 | ... | 0.002593 | 0 | 0 |
| 28 | 0.016887 | 0.016870 | 0.014460 | 0.012050 | ... | 0.002410 | 0 | 0 |
| 29 | 0.015667 | 0.015647 | 0.013411 | 0.011176 | ... | 0.002235 | 0 | 0 |
| ... | ... | ... | ... | ... | ... | ... | .. | .. |
| 8161 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8162 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8163 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8164 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8165 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8166 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8167 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8168 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8169 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8170 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8171 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8172 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8173 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8174 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8175 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8176 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8177 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8178 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8179 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8180 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |

| | | | | | | | | |
|------|----------|----------|----------|----------|-----|----------|---|---|
| 8181 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8182 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8183 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8184 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8185 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8186 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8187 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8188 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8189 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |
| 8190 | 0.000000 | 0.000000 | 0.000000 | 0.000000 | ... | 0.000000 | 0 | 0 |

| | 89 | 90 | 91 | 92 | 93 | 94 | 95 |
|------|----------|----|----|----|----------|----------|----------|
| 0 | 0.009020 | 0 | 0 | 0 | 0.027061 | 0.027061 | 0.027061 |
| 1 | 0.008833 | 0 | 0 | 0 | 0.026500 | 0.026500 | 0.026500 |
| 2 | 0.008633 | 0 | 0 | 0 | 0.025900 | 0.025900 | 0.025900 |
| 3 | 0.008421 | 0 | 0 | 0 | 0.025264 | 0.025264 | 0.025264 |
| 4 | 0.008199 | 0 | 0 | 0 | 0.024597 | 0.024597 | 0.024597 |
| 5 | 0.007967 | 0 | 0 | 0 | 0.023900 | 0.023900 | 0.023900 |
| 6 | 0.007726 | 0 | 0 | 0 | 0.023179 | 0.023179 | 0.023179 |
| 7 | 0.007478 | 0 | 0 | 0 | 0.022435 | 0.022435 | 0.022435 |
| 8 | 0.007225 | 0 | 0 | 0 | 0.021674 | 0.021674 | 0.021674 |
| 9 | 0.006966 | 0 | 0 | 0 | 0.020898 | 0.020898 | 0.020898 |
| 10 | 0.006703 | 0 | 0 | 0 | 0.020110 | 0.020110 | 0.020110 |
| 11 | 0.006438 | 0 | 0 | 0 | 0.019315 | 0.019315 | 0.019315 |
| 12 | 0.006172 | 0 | 0 | 0 | 0.018515 | 0.018515 | 0.018515 |
| 13 | 0.005905 | 0 | 0 | 0 | 0.017714 | 0.017714 | 0.017714 |
| 14 | 0.005638 | 0 | 0 | 0 | 0.016915 | 0.016915 | 0.016915 |
| 15 | 0.005373 | 0 | 0 | 0 | 0.016120 | 0.016120 | 0.016120 |
| 16 | 0.005111 | 0 | 0 | 0 | 0.015333 | 0.015333 | 0.015333 |
| 17 | 0.004852 | 0 | 0 | 0 | 0.014556 | 0.014556 | 0.014556 |
| 18 | 0.004597 | 0 | 0 | 0 | 0.013792 | 0.013792 | 0.013792 |
| 19 | 0.004348 | 0 | 0 | 0 | 0.013043 | 0.013043 | 0.013043 |
| 20 | 0.004103 | 0 | 0 | 0 | 0.012310 | 0.012310 | 0.012310 |
| 21 | 0.003865 | 0 | 0 | 0 | 0.011596 | 0.011596 | 0.011596 |
| 22 | 0.003634 | 0 | 0 | 0 | 0.010903 | 0.010903 | 0.010903 |
| 23 | 0.003410 | 0 | 0 | 0 | 0.010231 | 0.010231 | 0.010231 |
| 24 | 0.003194 | 0 | 0 | 0 | 0.009582 | 0.009582 | 0.009582 |
| 25 | 0.002985 | 0 | 0 | 0 | 0.008956 | 0.008956 | 0.008956 |
| 26 | 0.002785 | 0 | 0 | 0 | 0.008355 | 0.008355 | 0.008355 |
| 27 | 0.002593 | 0 | 0 | 0 | 0.007780 | 0.007780 | 0.007780 |
| 28 | 0.002410 | 0 | 0 | 0 | 0.007230 | 0.007230 | 0.007230 |
| 29 | 0.002235 | 0 | 0 | 0 | 0.006706 | 0.006706 | 0.006706 |
| ... | ... | .. | .. | .. | ... | ... | ... |
| 8161 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8162 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8163 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8164 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8165 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |

| | | | | | | | |
|------|----------|---|---|---|----------|----------|----------|
| 8166 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8167 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8168 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8169 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8170 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8171 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8172 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8173 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8174 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8175 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8176 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8177 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8178 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8179 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8180 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8181 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8182 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8183 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8184 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8185 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8186 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8187 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8188 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8189 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |
| 8190 | 0.000000 | 0 | 0 | 0 | 0.000000 | 0.000000 | 0.000000 |

[8191 rows x 94 columns]

We wish to find an expression

$$A_3 * \text{Spectrum}(Z = 3) + A_4 * \text{Spectrum}(Z = 4) + \dots + A_{95} * \text{Spectrum}(Z = 95) = \text{Spectrum}(Data) \quad (1)$$

such that the A_i represents how much of the $Z = i$ is present in the sample. We have an abundance of data, as each point represents a value that, as a scalar, must satisfy the above equation. Thus we convert this equation into its matrix form $\mathbf{A}\vec{x} = \vec{b}$. Each column of the above matrix (dataframe) represents the column of the coefficient matrix, \vec{x} is the column vector with elements A_i , and \vec{b} is the column vector describing the data obtained.

Because we have an overdetermined system, we do a linear least squares fitting to find the best 'fit'. However, we must add the constraint that $A_i \geq 0$, since we can't have spectral subtraction. For this, we use FORTRAN's non-negative least squares solver, implemented in `scipy.optimize.nnls`

```
In [42]: Z_dict = pd.read_fwf('X-Ray_database', names = ['Z', 'Element'], index_col = 0, skiprows=
# Z_dict['Element'].astype('str')
# print(Z_dict.iloc[16])
Z_dict.iloc[8]['Element']
Z_dict
```

Out[42]:

| | Element |
|-----|---------|
| Z | |
| 3 | Li |
| 4 | Be |
| 5 | B |
| 6 | C |
| 7 | N |
| 8 | O |
| 9 | F |
| 10 | Ne |
| 11 | Na |
| 12 | Mg |
| 13 | Al |
| 14 | Si |
| 15 | P |
| 16 | S |
| 17 | Cl |
| 18 | Ar |
| 19 | K |
| 20 | Ca |
| 21 | Sc |
| 22 | Ti |
| 23 | V |
| 24 | Cr |
| 25 | Mn |
| 26 | Fe |
| 27 | Co |
| 28 | Ni |
| 29 | Cu |
| 30 | Zn |
| 31 | Ga |
| 32 | Ge |
| ... | ... |
| 66 | Dy |
| 67 | Ho |
| 68 | Er |
| 69 | Tm |
| 70 | Yb |
| 71 | Lu |
| 72 | Hf |
| 73 | Ta |
| 74 | W |
| 75 | Re |
| 76 | Os |
| 77 | Ir |
| 78 | Pt |
| 79 | Au |
| 80 | Hg |

```

81      Tl
82      Pb
83      Bi
84      Po
85      At
86      Rn
87      Fr
88      Ra
89      Ac
90      Th
91      Pa
92      U
93      Np
94      Pu
95      Am

```

```
[93 rows x 1 columns]
```

```
In [43]: A = spectral_data[spectral_data.columns[1:]].as_matrix()
```

```


def composition(data):
    comps, resid = optimize.nnls(A,data)
    comps /= np.sum(comps)
    return comps, resid

def plot_composition(data, title, threshold):
    comps = composition(data)[0]
    important_idx = np.where(comps>threshold)[0]
    fig, ax = plt.subplots(figsize=(10,10))
    ax.bar(np.linspace(3,95, 93), comps*100, align = 'center')
    ax.set_xlabel('Atomic Number Z')
    ax.set_ylabel('\% Composition')
    for i in important_idx:
        ax.annotate('{:.1f}\% {}'.format(comps[i]*100, Z_dict.iloc[i]['Element']),\
                    xy = (i+3,comps[i]*100), xytext = (i+3,comps[i]*100+5),\
                    arrowprops=dict(facecolor='black', shrink=0.05) )
    ax.set_title(title)
    fig.tight_layout()


# Lets test this on the known spectrum for various elements as given by the website ref


plot_composition(spectral_data[8], 'Test Plot', .02)
# plot_composition(spectrograph('SumCrossXRF/Calibration/7Mar_Cu1.txt').get_counts(), 'Cu')
# plot_composition(spectrograph('SumCrossXRF/Data/7Mar_Gold.txt').get_counts(), 'Gold Fo')


```





```
In [59]: # Now lets look at the composition of all the data runs
for file in sorted(os.listdir('SumCrossXRF/Calibration/')):
    # print(file.partition('.')[0].partition('_')[2])
    spect = spectrograph('SumCrossXRF/Calibration/{}'.format(file))
    plot_composition(spect.get_counts(),file.partition('.')[0].partition('_')[2],.02)
```



















```
In [44]: # Now lets look at the composition of all the data runs
for file in sorted(os.listdir('SumCrossXRF/Data/')):
    # print(file.partition('.')[0].partition('_')[2])
    spect = spectrograph('SumCrossXRF/Data/{}'.format(file))
    plot_composition(spect.get_counts(),file.partition('.')[0].partition('_')[2],.02)
```
































Clearly, there are some populations in the composition that make me question the fit. For example, in the tin calibration samples, the dominant signal comes from sodium, which is more likely a false artefact from the data, as the sample seemed to be a reasonably well refined material.