

Laser Cooling and Trapping

February 27, 2018

```
In [2]: import numpy as np
import matplotlib.pyplot as plt
import pandas as pd
import scipy.optimize as optimize
from scipy.optimize import curve_fit
from scipy import signal

def log_growth(t, V_0, tau, t_0, V_offset):
    return V_0*(1-np.exp(-(t - t_0)/tau)) + V_offset
```

1 13 Feb 2017

SETUP

Going through equipment manuals and ensuring proper setup of equipment * Removed RF output to RF input cable between laser controller and laser servo. The equipment manual states that this is for peak lock only. * Connected Temp Servo input/output as per recommendation of manual for greater long term stability. * Tuned lasers using controller until trapping laser is on redshifted side of ^{87}Rb D_2 $5^2P_{3/2}$ $F' = 3$ line, and pump was on redshifted side of ^{87}Rb D_1 $5^2P_{3/2}$ $F' = 2$ line. * Locked lasers using the servo by adjusting the DC offset until the desired portion of each line was at 0 V, and then locked the lasers. * Confirmed the lasers are not hopping by ensuring the LD OUT line remained less than 50 mV.

2 15 Feb 2017

SETUP * Performed steps mentioned in previous section * We adjusted the mirrors on the table to ensure the beams all crossed in the center of the chamber. Because the beams are expanded, it is not as difficult to adjust, but we made sure the return beam reflected directly backwards as much as possible by checking for laser light around the perimeter of the primary director. * This Cold Quanta MOT has a different back panel than the manual. KEEP ION PUMP ON. We adjusted the dispenser current to approximately 3.3 mA in keeping with the Cold Quanta manual. However, there are no controls for the coils. * We found a power supply controlling the coils, and turned it on. We adjusted the current and voltage of the power supply for the coils to 1 V and 1 mA. * After turning on the dispenser for approximately 30 s, and toggling the polarization of the trapping beams, we noticed a glowing light on the CCTV observing the trap. * We blocked one of the trapping beams and noticed the light went away. After removing the block, the light came back over the course of several seconds. Jun suggested that this is a way to determine the stability of

our trap. * by adjusting the current and voltage of the coil power supply, we noticed that as we lowered these values, the ions flattened out perpendicular to the coil axis. * Over the next few days, we will devise several experiments and observations to make so we can have fun with this experiment * Is there a quantitative way to measure the brightness of the trap? Doing so could allow us to measure the filling time of the trap

3 20 Feb 2017

SETUP * We performed the same setup as before and immediately obtained a MOT. We will now build an imaging system to collect the light from the trapped atoms and measure the intensity of this light vs time. * Our MOT is ≈ 5.1 cm from the nearest we can place the lens. We use the below program to determine the optimal placement of lens and collector in the ≈ 220 mm available to us. * We found a 25 mm focal length lens, which should optimize the location of the image to the space constraints of our table.

```
In [3]: def lenseq(f,xo):
        """
        Given a lens focal length and location (giving xo), returns the image location xi
        """
        return xo*f/(xo-f)

        print(lenseq(25, 51))

49.03846153846154
```

- So our image should be approximately 5 cm from the lens.

4 22 Feb 2017

- It is difficult to obtain data from our photodetector, so we placed a camera in the place of the detector to optically inspect the what the photodetector is 'seeing'.
- We noticed that there is a great deal of light coming directly from the mirrors. These mirrors have a reasonably large angular separation from the MOT, so we need to place a tube around the detector to prevent these sources from entering the detector. We MacGyver'd a thin cardboard package and rolled it into a cylinder to achieve this blocking.
- Furthermore, we noticed that there is a significant amount of light from the laser scattering off the walls of the rubidium cell. We need something to limit the 'close-angle' light from entering the detector. We found an iris that we place close to the lens so we can limit the amount of scattering light that can enter the detector.
- It is difficult to 'aim' the photodetector to make sure it is 'seeing' the MOT. It took a great deal of adjustment to get the detector centered on the image.
 - In order to make sure we were seeing only the MOT, we would adjust the position with the MOT on, all while watching the signal from the detector.
 - Once we found a maximum, we would block all beams from entering the MOT area to ensure that the signal dropped, as the MOT should fail and so should any other laser light sources.


```

self.ax = plt.subplot(111)
for i in np.arange(len(plotdata_list)):
    plotdata_list[i].plot(kind = 'line', x = plotdata_list[i].columns[0], y = plotdata_list[i].columns[1], color = 'red', ax = self.ax, label = plotdata_list[i].columns[1], col
self.ax.set_xlim(self.data[self.data.columns[0]].min(), self.data[self.data.columns[1]].max())
self.ax.legend(loc = 'best')
self.fig.set_figheight(15)
self.fig.set_figwidth(15)
if set_lim:
    self.ax.set_ylim(self.data[self.data.columns[1]].min()-10, self.data[self.data.columns[1]].max()+10)
plt.show()

def zoom(self, dataset, lbound, rbound):
    asarray = np.array(dataset['t (s)'])
    l_idx = np.searchsorted(asarray, lbound)
    r_idx = np.searchsorted(asarray, rbound)

    return dataset[l_idx:r_idx]

def filter_data_butter(self, butterworth_order, cutoff_freq, pad):
    filtered_data = self.data.copy()
    b, a = signal.butter(butterworth_order, cutoff_freq)
    filtered_data[filtered_data.columns[1]] = signal.filtfilt(b,a, filtered_data[filtered_data.columns[1]],
                                                             padlen = pad)

    return filtered_data

def filter_data_gust(self, ellip_order, max_ripple, min_attenuation, cutoff_freq):
    filtered_data = self.data.copy()
    b, a = signal.ellip(ellip_order, max_ripple, min_attenuation, cutoff_freq)
    filtered_data[filtered_data.columns[1]] = signal.filtfilt(b,a, filtered_data[filtered_data.columns[1]],
                                                             method = "gust")

    return filtered_data

def logistic_fit(self, data, initial_guess):
    """
    Fits data to the curve  $V(t) = V_0 * (1 - \exp(-(t-t_0)/\tau)) + V_{offset}$ 
    Requires intial guess to be of the form [V_0, tau, t_0, V_offset]
    data must have two columns: first column is time, second in voltage
    """
    best_vals, covar = curve_fit(log_growth, data[data.columns[0]], data[data.columns[1]])
    dTau = best_vals[1]
    errTau = np.sqrt(np.diag(covar))[1]

    return best_vals, covar

def make_bestfit(self, function, x_data, fit_vals, val_err):
    report_Tau, report_Err = self.reported_values(fit_vals[1], np.sqrt(np.diag(val_err)))

```

```

xmin = self.data[self.data.columns[0]].min()
xmax = self.data[self.data.columns[0]].max()
x_data = np.linspace(xmin, xmax, 1000)
#     bestfit_data = pd.DataFrame([x_data] , [function(x_data, *fit_vals)],\
#     columns = [self.data.columns[0], 'Best Fit (mV)'])
bestfit_data = pd.DataFrame({self.data.columns[0]:x_data})
bestfit_data[ 'Best Fit (mV)\n$ \\tau $ = {}$\\pm${}$'.format(report_Tau, report_E

return bestfit_data

def reported_values(self, value, error):
    a = int(np.floor(np.log10(np.abs(error))))
    report_Err = round(error, -a)
    report_Val = round(value, -a)


    return report_Val, report_Err

def best_w_data(self, data_list, best_guess_list, color_list = ['r','b','r', 'g', 'r']
    """
    Given a list of data and the associated best guesses for the fit,\
    plots the curves and best fits for those plots
    """
    data_w_best = []
#     data_w_best = [test.data]
    for i in np.arange(len(data_list)):
        data_w_best.append(data_list[i])
        if i != 0:
            best_vals, val_err = self.logistic_fit(data_list[i], best_guess_list[i])
            data_w_best.append(self.make_bestfit(log_growth, data_list[i][data_list[


    self.plot_data(data_w_best,color_list, set_lim=True)

In [5]: test0 = MOT_data('./Data/22FEB2018/tek0000CH3.csv', 't (s)', 'Photodetector Voltage (mV)')
# test0.plot_data([test0.data])
test0.plot_data([test0.zoom(test0.data,-75,100)])


```




```
In [6]: test0.best_w_data([test0.data, test0.zoom(test0.data, -78,100)] ,\
                        [[],[3, 75., -80.,-20]])
```




```
In [7]: test0.best_w_data([test0.filter_data_gust(5,0.001,120,.001),test0.zoom(test0.filter_data
      [],[3, 75., -80.,-20])])
```




```
In [8]: test1 = MOT_data('./Data/22FEB2018/tek0001CH3.csv', 't (s)', 'Photodetector Voltage (mV)')
# test1.plot_data([test1.data])
test1.plot_data([test1.zoom(test1.data,-75,100)])
```


```
In [9]: test1.best_w_data([test1.data, test1.zoom(test1.data, -78,100)] ,\n                        [[],[3, 75., -80.,-20]])
```




```
In [10]: test1.best_w_data([test1.filter_data_gust(5,0.001,120,.001),test1.zoom(test1.filter_data_gust(5,0.001,120,.001),[3, 75., -80.,-20])])
```




```
In [11]: test3 = MOT_data('./Data/22FEB2018/tek0003CH3.csv', 't (s)', 'Photodetector Voltage (mV)')
# test3.plot_data([test3.data])
test3.plot_data([test3.zoom(test0.data,-72,100)])
```




```
In [12]: test3.best_w_data([test3.data, test3.zoom(test3.data, -72,100)] ,\
                           [[],[3, 75., -70.,-19]])
```




```
In [13]: test3.best_w_data([test3.filter_data_gust(5,0.001,120,.001),test1.zoom(test3.filter_data_gust(5,0.001,120,.001),[3, 75., -70.,-19])])
```




```
In [14]: test4 = MOT_data('./Data/22FEB2018/tek0004CH3.csv', 't (s)', 'Photodetector Voltage (mV)')
# test4.plot_data([test4.data])
test4.plot_data([test4.zoom(test4.data,-79,20),test4.zoom(test4.data,25,100)])
```




```
In [15]: test4.best_w_data([test4.data, test4.zoom(test4.data,-79,20), test4.zoom(test4.data,25,
        [ [],[3, 35., -80.,-19], [3, 35., 20.,-19]])
```




```
In [16]: test4.best_w_data([test4.filter_data_gust(5,0.001,120,.001),test4.zoom(test4.filter_data_gust(5,0.001,120,.001),25,96)] ,\
                           test4.zoom(test4.filter_data_gust(5,0.001,120,.001),25,96)) ,\
                           [[],[3, 35., -80.,-19], [3, 35., 20.,-19]])
```


```
In [17]: test5 = MOT_data('./Data/22FEB2018/tek0005CH3.csv', 't (s)', 'Photodetector Voltage (mV)')
test5.plot_data([test5.zoom(test5.data, -79,22), test5.zoom(test5.data,22.35,97.8)])
```




```
In [18]: test5.best_w_data([test5.data, test5.zoom(test5.data, -79,22), test5.zoom(test5.data,22,100,20)],  
                           [[], [10, 50., -80., -20], [10, 50., 20., -20]])
```




```
In [19]: test5.best_w_data([test5.filter_data_gust(5,0.001,120,.001),test5.zoom(test5.filter_data_gust(5,0.001,120,.001),23,96)] ,\
                             test5.zoom(test5.filter_data_gust(5,0.001,120,.001),23,96)) ,\
                             [[],[10, 50., -80.,-20], [10, 50., 20.,-20]])
```




```
In [20]: # Compare fits of the different filters on each peak
test5.best_w_data([test5.zoom(test5.filter_data_gust(5,0.001,120,.002), -79,22),\
test5.zoom(test5.filter_data_gust(5,0.001,120,.001), -79,22),\
test5.zoom(test5.filter_data_butter(10,.025, 150), -79,22)],\
[[10, 50., -80.,-20], [10, 50., -80.,-20],[10, 50., -80.,-20]])
```




```
In [21]: test5.plot_data([test5.filter_data_butter(10,.025, 150), test5.filter_data_gust(5,0.001
```




```
In [22]: test6 = MOT_data('./Data/22FEB2018/tek0006CH3.csv', 't (s)', 'Photodetector Voltage (mV)')
# test6.plot_data([test6.data])#test6.zoom(test6.data,22.35,97.8)
test6.plot_data([test6.zoom(test6.data, -78,20), test6.zoom(test6.data,23,97.8)])
```




```
In [23]: test6.best_w_data([test6.data, test6.zoom(test6.data,-78,20), test6.zoom(test6.data,21,
[[],[5, 50., -78.,-20], [5, 50., 21.,-20]])
```




```
In [24]: test6.best_w_data([test6.filter_data_gust(5,0.001,120,.001),test6.zoom(test6.filter_data_gust(5,0.001,120,.001),21,95)],\
                             test6.zoom(test6.filter_data_gust(5,0.001,120,.001),21,95)) ,\
                             [[],[5, 50., -78.,-20], [5, 50., 21.,-20]])
```


```
In [25]: test7 = MOT_data('./Data/22FEB2018/tek0007CH3.csv', 't (s)', 'Photodetector Voltage (mV)')
# test7.plot_data([test7.data])#test6.zoom(test6.data,22.35,97.8)
test7.plot_data([test7.zoom(test7.data, -72,25), test7.zoom(test7.data,27,92)])
```




```
In [26]: test7.best_w_data([test7.data, test7.zoom(test7.data,-72,25), test7.zoom(test7.data,27,
[[],[5, 50., -78.,-20], [5, 50., 21.,-20]])
```




```
In [27]: test7.best_w_data([test7.filter_data_gust(5,0.001,120,.001),test7.zoom(test7.filter_data_gust(5,0.001,120,.001),27,92)] ,\
                           test7.zoom(test7.filter_data_gust(5,0.001,120,.001),27,92)] ,\
                           [[],[5, 50., -78.,-20], [5, 50., 21.,-20]])
```



```
In [28]: test8 = MOT_data('./Data/22FEB2018/tek0008CH3.csv', 't (s)', 'Photodetector Voltage (mV)')
# test8.plot_data([test8.data])
test8.plot_data([test8.zoom(test8.data, -76,20), test8.zoom(test8.data,27,98)])
```



```
In [29]: test8.best_w_data([test8.data, test8.zoom(test8.data, -76,20), test8.zoom(test8.data,27
          [, [5, 50., -78., -20], [5, 50., 21., -20]])
```



```
In [30]: test8.best_w_data([test8.filter_data_gust(5,0.001,120,.001),test8.zoom(test8.filter_data_gust(5,0.001,120,.001),27,97)],\
                             test8.zoom(test8.filter_data_gust(5,0.001,120,.001),27,97)) ,\
                             [[],[5, 50., -78.,-20], [5, 50., 21.,-20]])
```

